

## YAD SARAH Powered Mobility Aids Approval

EL AL Dangerous good team approve the carriage of the following scooter/wheelchair on El Al flights:

Mobility aid Model	Lithium Battery Model	Wh /(V)X(Ah)	Manufacture Battery
52016/52092	FY24-6 or FDBO1	144	CHANGZHOU FENGYU NEW ENERGY TECH.
52013/52134	N3458 (18650-7S4P)	240	HENGDIAN GROUP DMEGC MAGNETICS
52134	LI18650 (7S3P)	124.8	HENGDIAN GROUP DMEGC MAGNETICS
52134	FK3458-10D	240	YONGKANG FOTEC ELECTRONIC TECHNOLOGY.

### The battery designed to be removed by the user.

Carriage requirements will be verified at the time of Check-in:

1. The battery must be removed from the scooter by the passenger.
  2. The battery must be protected by placing the battery in a separate bag or by taping over exposed terminals.
  3. The battery must be carried in the passenger cabin by the passenger.
- If any of the above details have been changed, you must inform El Al.
  - In order for EL AL out base stations to prepare for your boarding, it is very important to update EL AL at least 48 hours before your check in. please send this approval with the proper model number to [MEDADEP@ELAL.CO.IL](mailto:MEDADEP@ELAL.CO.IL)
  - more information about traveling with mobility aids and dangerous goods can be found on El Al website: [Information on batteries and dangerous substances | EL AL](#)