

YAD SARAH Non-Spillable Powered Mobility Aid Approval

EL AL Dangerous good team approve the carriage of the following scooter on El Al flights:

Mobility aid Model	Non Spillable – VRLA Gel Battery Model	Wh /(V)X(Ah)	Manufacture Battery
52001	6-DZF-20	12VX12Ah	CHILWEE



Carriage requirements will be verified at the time of Check-in:

1. The battery terminals must be protected from short circuits. (Enclosed within a battery container).
 2. The battery must be securely attached to the scooter and the electrical circuits are isolated following the manufacturer's instructions.
 3. The scooter with the battery must be carried in the cargo compartment.
 4. The passenger can carry maximum of one spare battery.
- If any of the above details have been changed, you must inform El Al.
 - In order for EL AL out base stations to prepare for your boarding, it is very important to update EL AL at least 48 hours before your check in. please send this approval with the proper model number to MEDADEP@ELAL.CO.IL
 - more information about traveling with mobility aids and dangerous goods can be found on El Al website: [Information on batteries and dangerous substances | EL AL](#)